



Product Information

TRM Staysound

Name:	TRM Staysound	
Pack Sizes/Codes:	1.5kg (TRMSSP1.5) 5kg (TRMSSP5) 11.35kg (TRMSSP11.35)	
Main Points:	<ul style="list-style-type: none"> • Kaolin clay to soothe and cool tired legs • Ready to use, easy to apply and long lasting formula 	
Status:	Vet OTC	
Ingredients:	<ul style="list-style-type: none"> • Montmorillonites • Aluminium Silicate • Sodium Borate • Glycerine • Natural Thyme Oil • Purified water 	
Indications:		
Action:		
Administration:	<p>Routine use: apply directly to the knees, fetlocks, shins and tendons after periods of training. Spread about 6mm thick</p> <p>Prolonged effect: Wrap with moistened brown paper and bandage. Repeat as necessary after 24 hours</p> <p>Remove by hosing with water.</p>	
Dosage:		
Special Precautions:		
Storage:	N/A	
Registration:	Exempt	
Category:	Mobility and joint care	



Animal Health Direct Limited

PO Box 8015 | Havelock North | 4157 | t: 06 873 3611

e: sales@ahdltd.co.nz | ahdltd.co.nz

